NEOVALVE®

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NEOVALVE® Sterile Suction Valve

Please read this leaflet carefully before using this medical device because it contains important information.
If you have more questions, please contact Asept InMed.

1. DESCRIPTION

The range is composed of 2 types of sterile suction control valve with or without suction tube:

- Sterile suction control valve adult, with or without suction tube with different lengths;
- Sterile suction control valve pediatric, with or without suction probe with different lengths.

2. INDICATIONS

This device is intended for aspiration of tracheal and bronchial mucus in order to clean the airway and facilitate breathing.
Use of a probe connected to the suction valve is required for performing this suction.

3. CONTRAINDICATIONS

Suction through the nose is prohibited.

4. COMPLICATIONS

On a conscious person, it can lead to risks of vomiting (nausea reflex).
Respiratory risks (desaturation or even hypoxia, asphyxia of the patient).
Risks of bleeding from the probe connected to the suction valve are possible.

5. WARNINGS

Important / very thick secretions can interfere with proper functioning of the valve.
Do not use if the product is damaged, present a defect or is not functional.
Reuse or cleaning the device can alter their structural and mechanical properties. It can lead to infection or other illness/injury.
Poly Medicure will not be responsible for any direct incidental or consequential damages resulting from reuse of devices.
Use according indications.

6. USE PRECAUTIONS

Device must be used exclusively by health professional.
The aspiration should be done as clean as possible and over a short period (10 seconds for adults, 4 seconds for pediatric) in order to limit risks of traumatic injuries and infections.

7. HYGIENE PRECAUTIONS

Respect the general rules of hygiene and asepsis of health establishments.
For user safety, it could be necessary to use sterile gloves, glasses, mask during suction process.

8. FREQUENCY OF USE

Single patient device.
Suction should be done every 4 to 8 hours or at every service change.
The device should be changed every 24h to 48 hours.

9. MODALITY OF USE

It is recommended to wash or disinfect your hands before and after use.

- 1- Connect the device on the hook side to the suction probe (except for suction tube already connected) complying the direction of assembly. Firmly connect the components during assembly (except for suction tube already connected) ;
- 2- Perform a functional test (adjusting the vacuum at -500mb):
 - o Check the absence of air leak on the suction side ;
 - o Close the hole of the cover with finger and check the correct suction ;
- 3- In case of air leakage or whistling sound, check the integrity of the probes and connections and increase the vacuum until a correct process (check the seal of the container, cap improperly closed, device improperly connected) ;
- 4- In case of air leakage or improper suction, change the system and repeat the test ;
- 5- To limit any lack of oxygen, each suction should exceed 10 seconds for adults and 4 seconds for children and new born. Change the device in case of important/abundant secretion.

Device should be rince with sterile water or with a rincing solution after use. Container should be replace after 24h of use.
Suction probe should be change after use.

10. STORAGE AND TRANSPORT CONDITIONS

Must be conserved, stored and transported avoiding excessive heat, protecting from direct sunlight and moisture.

11. DISPOSABLE CONDITIONS

Used devices must be eliminated according to the hygiene protocol of the health facility and to the actual regulation.

12. PERFORMANCES

Performances of the device are technical.




13. MATERIAL

Valves and pediatric connectors are in Styren Butadiene copolymer and the membrane is in rubber Nitrile.
The suction tube is in Polyvinylchloride.

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These devices are unitary packaged and sterilized with ethylene oxyd..
Do not use if the primary packaging is opened or damaged.
Expiring delay: 5 years
Do not use if expired.

15. SYMBOLS

| | | | |
|---|------------------------------|---|------------------------------------|
|  | Manufacturer |  | Single use |
|  | Manufacturing date |  | Do not use if packaging is damaged |
|  | Warning |  | Reference |
|  | Sterilized with ethylen oxyd |  | Batch number |
|  | Keep dry |  | Expiry date |
|  | Keep away from sunlight |  | European Representative |

16. Administrative information

Date of modification of the leaflet: April 2021.

If you have any further questions or need information, please contact us at the following address:

 Poly Medicure
Plot No. 115-116, Sector 59, HSIDC Indi. Area,
Faridabad-121004, (Haryana) INDIA
www.polymedicure.com
Customer Care No. : +91-129-4287053

Imported by:
ASEPT INMED
9, Avenue Mercure 31130 Quint Fonsegrives France
Tél : +33(0)5 62 57 69 00 Fax : +33(0)5 62 57 69 01
www.aseptinmed.fr

 Obelis s.a
Boulevard Général Wahis 53
1030 Brussels, BELGIUM
Tel : + (32) 2.732.59.54
E-mail : mail@obelis.net