# microlife®

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## Références

- [1] Nichols M TN, Luengo-Fernandez R, Leal J, Gray A, Scarborough P, Rayner M European Cardiovascular Disease Statistics 2012. European Heart Network, Brussels, European Society of Cardiology, Sophia Antipolis 2012.
- [2] Boriani G, Laroche C, Diemberger I, Fantecchi E, Popescu MI, Rasmussen LH et al.: Asymptomatic atrial fibrillation: clinical correlates, management, and outcomes in the EORP-AF Pilot General Registry. Am J Med 2015; 128:509-518 e502.
- [3] Zoni-Berisso M, Lercari F, Carazza T, Domenicucci S: Epidemiology of atrial fibrillation: European perspective. Clin Epidemiol 2014; 6:213-220.
- [4] Sanna T, Diener HC, Passman RS, Di Lazzaro V, Bernstein RA, Morillo CA et al.: Cryptogenic stroke and underlying atrial fibrillation. N Engl J Med 2014; 370:2478-2486.
- [5] Gladstone DJ, Spring M, Dorian P, Panzov V, Thorpe KE, Hall J et al.: Atrial fibrillation in patients with cryptogenic stroke. N Engl J Med 2014; 370:2467-2477.
- [6] Camm AJ, Lip GY, De Caterina R, Savelieva I, Atar D, Hohnloser SH et al.: 2012 focused update of the ESC Guidelines for the management of atrial fibrillation: an update of the 2010 ESC Guidelines for the management of atrial fibrillation--developed with the special contribution of the European Heart Rhythm Association. Europace 2012; 14:1385-1413.
- [7] Hart RG, Benavente O, McBride R, Pearce LA: Antithrombotic therapy to prevent stroke in patients with atrial fibrillation: a meta-analysis. Ann Intern Med 1999; 131:492-501.
- [8] Ruff CT, Giugliano RP, Braunwald E, Hoffman EB, Deenadayalu N, Ezekowitz MD et al.: Comparison of the efficacy and safety of new oral anticoagulants with warfarin in patients with atrial fibrillation: a meta-analysis of randomised trials. Lancet 2014; 383:955-962.
- [9] Verberk WJ, Omboni S, Kollias A, Stergiou GS: Screening for atrial fibrillation with automated blood pressure measurement: Research evidence and practice recommendations. Int J Cardiol 2016; 203:465-473.
- [10] Wiesel J, Wiesel D, Suri R, Messineo FC: The use of a modified sphygmomanometer to detect atrial fibrillation in outpatients. Pacing Clin Electrophysiol 2004; 27:639-643.
- [11] Stergiou GS, Karpettas N, Protopogrou A, Nasothimiou EG, Kyriakidis M: Diagnostic accuracy of a home blood pressure monitor to detect atrial fibrillation. J Hum Hypertens 2009; 23:654-658.
- [12] Wiesel J, Fitzig L, Herschman Y, Messineo FC: Detection of atrial fibrillation using a modified microlife blood pressure monitor. Am J Hypertens 2009; 22:848-852.
- [13] Wiesel J, Abraham S, Messineo FC: Screening for asymptomatic atrial fibrillation while monitoring the blood pressure at home: trial of regular versus irregular pulse for prevention of stroke (TRIPPS 2.0). Am J Cardiol 2013; 111:1598-1601.
- [14] Kearley K, Selwood M, Van den Bruel A, Thompson M, Mant D, Hobbs FR et al.: Triage tests for identifying atrial fibrillation in primary care: a diagnostic accuracy study comparing single-lead ECG and modified BP monitors. BMJ Open 2014; 4:e004565.
- [15] Wiesel J, Arbesfeld B, Schechter D: Comparison of the Microlife blood pressure monitor with the Omron blood pressure monitor for detecting atrial fibrillation. Am J Cardiol 2014; 114:1046-1048.
- [16] Gandolfo C, Balestrino M, Bruno C, Finocchi C, Reale N: Validation of a simple method for atrial fibrillation screening in patients with stroke. Neurol Sci 2015; 36:1675-1678.
- [17] Chan PH, Wong CK, Pun L, Wong YF, Wong MM, Chu DW et al.: Head-to-Head Comparison of the AliveCor Heart Monitor and Microlife WatchBP Office AFIB for Atrial Fibrillation Screening in a Primary Care Setting. Circulation 2017; 135:110-112.
- [18] Chan PH, Wong CK, Pun L, Wong YF, Wong MM, Chu DW et al.: Diagnostic performance of an automatic blood pressure measurement device, Microlife WatchBP Home A, for atrial fibrillation screening in a real-world primary care setting. BMJ Open 2017; 7:e013685.
- [19] Alpert BS, Quinn D, Gallick D: Oscillometric blood pressure: a review for clinicians. J Am Soc Hypertens 2014; 8:930-938.
- [20] NICE: WatchBP Home A for opportunistically detecting atrial fibrillation during diagnosis and monitoring of hypertension <http://guidance.nice.org.uk/MTG13>. 2013; Assessed 18 Aug. 2015.