



## NON-STERILE DELICATE STRIP PLASTER (BOY & GIRL) WITH PAD TREATED WITH ANTIBACTERIAL AGENT TECHNICAL INFORMATION

### **Definition**

- Delicate strip plaster for boys – in packs.

### **General characteristics**

- Ready-to-use strip plaster, consisting of a nonwoven fabric backing coated with hypoallergenic adhesive with high tolerability for skin and with a central absorbent pad treated with antibacterial agent (0.2% chlorhexidine digluconate) that controls the pad's micro-environment but does not come into contact with the wound.
- It has an extremely conformable adhesive backing that follows the contours of the skin to ensure patient comfort during application.
- The cross-direction extensibility of the nonwoven fabric of the backing makes the dressing suitable for application on joints and moving parts.
- The backing has a good degree of adhesiveness that ensures constant grip in time and a good resistance to changes in temperature.
- Very delicate removal and leaves no adhesive residues on skin.
- The perforated backing and the application of the adhesive in rows allow the proper flow of air and water vapour thus ensuring skin transpiration.
- The central pad is highly absorbent to protect any type of wound and absorb the various levels of exudate. It contains an antiseptic (0.2% chlorhexidine digluconate) that regulates the pad's micro-environment but does not come into contact with the wound.
- The part of the pad that comes into contact with the wound is constituted in a non-adherent net that reduces adhesion to the wound to a minimum.
- The pad and adhesive backing are covered with silicone-coated paper that can be easily removed before use.
- It is radiotransparent, which means that it is invisible to X-rays.
- Available in medium sizes and it is packed in a box containing 24 pieces (see assortment).



### **Indications**

- For the dressing of wounds in hospitals and outpatient clinics.
- In all situations that require dressing the wound with a fast technique not requiring secondary dressings.
- Recommended for dressing wounds on all types of skin and particularly for sensitive skin.
- Ideal for dressings of joints, parts involved in frequent movement or irregular surfaces of the body.

**Medical device according to Directive 93/42/EEC,**

**modified by Directive 2007/47/EC**

**Implementation with Legislative Decree no. 37 of 25-01-2010**

**CE 0373**

### **Class IIb**

- Intended use: non-sterile non-invasive surgical dressing device, with absorbent pad treated with antiseptic (0.2% chlorhexidine digluconate) that regulates the micro-environment of the pad in contact with the wound and with adhesive backing in contact with intact skin.

Pikdare srl considers the information herein as confidential and not disclosable. Pikdare reserves the exclusive right to modify the document according to the state-of-the-art and regulatory changes concerning the product.



## NON-STERILE DELICATE STRIP PLASTER (BOY & GIRL) WITH PAD TREATED WITH ANTIBACTERIAL AGENT TECHNICAL INFORMATION

### Compliances

- F.U./F.E. Vigenti - DIRETTIVA CEE 93/42
- ISO 2859 Piani di campionamento
- ISO 10993 Biocompatibilità
- Sistema Assicurazione qualità: certificato secondo UNI EN ISO 9001:2008
- Sistema Assicurazione qualità: certificato secondo EN ISO 13485:2012

### Technical specifications

SPECIFICATIONS	METHODS	VALUES
<b>ADHESIVE BACKING</b>		
Composition of adhesive backing		100% polyester nonwoven fabric
Colour of adhesive backing		Printed with coloured drawings
Weight of adhesive backing	EDANA 40-3/90	50 ±5 g/m <sup>2</sup>
Thickness of adhesive backing	EDANA 30-5/99	0,46 - 0,60 mm
Machine-direction breaking load (MD)	FU IX ed.	≥ 18N/cm
Machine-direction elongation at break (MD)	FU IX ed.	≤50%
Cross-direction breaking load (CD)	FU IX ed.	≥3.6 N/cm
Cross-direction elongation at break (CD)	FU IX ed.	≤170%
Machine-direction extensibility (MD)	FU IX ed.	Not extensible
Cross-direction extensibility (CD)	FU IX ed.	Extensible
Composition of adhesive mass		Solvent-based synthetic acrylic adhesive
Weight of adhesive mass	M.I.	40 ±2 g/m <sup>2</sup>
Weight of adhesive backing	M.I.	90 ±7 g/m <sup>2</sup>
Coating of adhesive mass		In rows
Adhesiveness on steel	FU IX ed.	5,0-10,0 N/cm
Moisture vapour transmission rate (MVTR)	FU IX ed.	>3.500 g/m <sup>2</sup> /24h/37°C
<b>ABSORBENT PAD</b>		
Composition of absorbent part		Non-woven fabric composed of 85% viscose and 15% polypropylene, pad treated with 0.2% chlorhexidine digluconate
Colour of the absorbent part		White
Composition of non-adherent net		Polyethylene film
Colour of non-adherent net		White
Total weight of pad	M.I.	106 ±10 g/m <sup>2</sup>
Absorbency	UNI EN 13726-1	6.8±2 g/100 cm <sup>2</sup>
<b>PROTECTIVE PAPER</b>		
Composition of protective paper		Silicone-coated paper
Colour of protective paper		White
Weight of protective paper	M.I.	62 ±3 g/m <sup>2</sup>

(MD: Machine Direction; CD: Cross Direction)

\* According to UNI EN 13726, the moisture vapour transmission rate and conformability of a primary dressing are the result of the combination of the layers composing it. These tests are carried out on the complete dressing: absorbent pad coupled with the adhesive backing.

Pikdare srl considers the information herein as confidential and not disclosable. Pikdare reserves the exclusive right to modify the document according to the state-of-the-art and regulatory changes concerning the product.



## NON-STERILE DELICATE STRIP PLASTER (BOY & GIRL) WITH PAD TREATED WITH ANTIBACTERIAL AGENT TECHNICAL INFORMATION

<b>Sterile</b>	<b>No</b>
<b>Latex content in the product</b>	<b>No</b>
<b>Latex content in the package</b>	<b>Yes, in the seal of the packaging paper</b>
<b>Validity</b>	<b>5 years in proper storage conditions</b>
<b>Storage and shelf life</b>	<b>Store in a cool, dry place</b>

### Possibility to sterilise the product

- Sterilisable, in appropriate packaging, by irradiation for a maximum bacterial load of <100 ufc/g.

### Number of possible resterilisations

/

### Biocompatibility

- The product has passed the biocompatibility tests required by standard ISO 10993: cytotoxicity, allergic sensitisation, skin irritation.
- Compatible with medicines and disinfectants normally used.

### Disposal

- In compliance with applicable laws.

### Warnings

- Ensure clean and dry hands before using the product.
- Store in a cool, dry place.
- Apply the adhesive backing on clean, intact and dry skin in order not to alter the adhesiveness of the adhesive backing.

### Quality testing

- Quality testing is carried out on the raw material, semi-finished product and finished product. Tests are carried out to assess conformity to Pharmacopoeia or internal reference specifications.

### Packaging

The labels and/or packaging bear the following information:

Product name
Qualitative and quantitative composition
Dimensions
Instructions for use and warnings
Manufacturer
Date and batch of manufacture
Validity
Storage
Item barcode

### Packaging material

- Self-sealing paper pouch (contains natural latex).
- Printed cardboard pack.
- The packaging and materials and inks are certified as non-toxic.

Pikdare srl considers the information herein as confidential and not disclosable. Pikdare reserves the exclusive right to modify the document according to the state-of-the-art and regulatory changes concerning the product.



## NON-STERILE DELICATE STRIP PLASTER (BOY & GIRL) WITH PAD TREATED WITH ANTIBACTERIAL AGENT TECHNICAL INFORMATION

### Manufacturer

➤ Pkdare S.r.l. - Via Saldarini Catelli 10 - 22070 Casnate con Bernate (Como) - Italy.

### Assortment

CODE	DIMENSIONS mm x mm	PACK	UNIT	CARTON
0202409000000	19X72	24 pcs	12 sc	192 sc
0202409100000	19X72	24 pcs	12 sc	192 sc

### Other useful information

➤ Additional information and operating procedures are contained in the technical dossiers filed with the PIKDARE Technical/Quality Assurance Departments.

**Technical Director**

PIKDARE S.R.L.  
Dr. MAURO CASSANI  
Medical Devices Department  
Technical Director

Pikdare srl considers the information herein as confidential and not disclosable. Pikdare reserves the exclusive right to modify the document according to the state-of-the-art and regulatory changes concerning the product.