

Specifications

Basic UDI-DI: 897564000733MH

1. Variants:

Item no.	Variant	GTIN / GS1 UDI-DI
2600	Omni Tray	00897564000351

2. Bill of materials:

ITEM NO.	PART NUMBER	REV	DESCRIPTION	MATERIAL	QTY
1	22-0061	A	NUT, HEX FLANGED, M6 X 1.0	CRS	2
2	22-0084	A	WASHER, M10, 22MM OD	CRS	1
3	22-0086	C	NUT, M10 X 1.5MM, NYLOCK	CRS	1
4	1000-013	A	NUT, LOW PROFILE, M8-1.25	PLAIN CARBON STEEL	2
5	1000-036	B	NUT, NYLOCK, M6, BLACK NICKEL	CRS	1
6	1000-069	A	BOLT, BH, M6 X 1.0 X 40MM, 10.9, BLACK ZINC	CRS	2
7	1000-070	A	WASHER, FLAT, M6, BLACK ZINC	CRS	1
8	1000-119	A	SCREW, WOOD	CRS	12
9	1000-246	A	BOLT, CARRIAGE, M6 X 1.0 X 40MM	CRS	1
10	2600-037	B	BAMBOO TRAY	BAMBOO	1
11	1000-270	A	SCREW, SET, M10 X 1.5 X 12MM	CRS	3
12	1000-275	A	WELDNUT, M10, WITH NOTCH	CRS	3
13	1000-290	1	BOLT, CARRIAGE, M6 X 1.0 X 15, WIDE HEAD	CRS	1
14	1000-291	A	SPRING PIN FOR 8.5MM HOLE	SPRING STEEL	1
15	1000-300	A	WASHER, SPLIT, M8, BLACK ZINC	CRS	2
16	1000-301	B	BOLT, BH, M8 X 1.25 X 45MM, GRADE 10.9	CRS	2
17	1000-5020	B	BOLT, UNION, FEMALE, M6 X 1.0 X 15MM	CRS	1
18	2600-004	D	EXTENDABLE UPRIGHT TUBE, 25MM	Q195	1
19	2600-006	D	HANDLE TUBE POST	SPCC Q195	1
20	2600-007	C	BASE LONG PLASTIC	ABS	1

21	2600-008	B	TOP PLASTIC	ABS	1
22	2600-010	D	BASE TUBE, 30MM	SPCC Q195	2
23	2600-011	A	CC BASE PLATE	SPCC Q195	2
24	2600-012	B	CC FOOT PAD	TPR	4
25	2600-013	A	CC RUBBER FOOT	TPR	4
26	2600-021	A	CAP, SQUARE END 30MM	PP	2
27	2600-022	B	CAP, SQUARE OPEN END 30MM	PP	3
28	2600-023	D	BASE EXTENSION TUBE	SPCC Q195	1
29	2600-024	B	CAP, INTERNAL SQUARE END 25MM	PP	3
30	2600-025	1	WELDMENT, GUSSET	SPCC Q195	4
31	2600-030	F	35MM LOWER UPRIGHT BRACKET	SPCC Q195	1
32	2600-031	H	30MM UPRIGHT TUBE	SPCC Q195	1
33	2600-032	C	BUSHING, TRAY ARM SWIVEL	POM	1
34	2600-033	D	SLEEVE, TRAY ARM PIVOT	SPCC Q195	1
35	2600-034	D	TRAY ARM TUBE	SPCC Q195	1
36	2600-036	B	CAP, SQUARE END, 25MM	PP	1
37	2600-038	E	TRAY BRACKET	CRS	1
38	2600-039	E	TABLE BRACKET PIN	CRS	1
39	2600-041	D	TRAY PIVOT	A380 / ADC12	1
40	2600-042	C	WASHER, TRAY PIVOT	POM	1
41	2600-043	E	BUSHING, TABLE SWIVEL	PP	1
42	2600-045	B	ORGANIZER, TRAY	PP	1
43	2600-046	B	CUP HOLDER, TRAY	PP	1
44	2600-047	B	CAP, DOMED END	PP	1
45	2600-048	B	BIN, UPRIGHT	PP	1
46	2600-049	C	HEIGHT ADJUSTABLE BAND	TPR	1

3. Specifications:

Variant	Product Weight	Weight Capacity of Handle	Weight Capacity of Tray	Handle Size	Tray Size	Height Adjustment of Handle
2600 Omni Tray	9 KG	113.4 KG	13.6 KG	11.4 CM X 12.7 CM	50.8 CM X 38 CM	72.4 CM- 92 CM