

Ref.: 80111S
 Fecha/date: 28-06-2023



E0015114 (EAN010C) (SP)
 LPP 6223534 (FR)



Descripción Description

1) Andador de 4 ruedas con asiento y cesta que permite el apoyo a personas con dificultades para caminar y así poder sostener y mantener el equilibrio mientras caminan.

2) Homologado producto sanitario en base a los requisitos generales de seguridad, funcionamiento y disposiciones del Reglamento (EU) 2017/745 del Parlamento Europeo y la norma específica: UNE-EN ISO 11199-2:2021.

3) Freno de mano con posibilidad de bloqueo de ruedas traseras.

4) Ruedas grandes con bloqueo 0-90-360°.

5) Empuñadura ajustable en altura.

6) Bolsa extraíble con capacidad para 25 litros.

7) Amplio soporte lateral ergonómico.

8) Plegado compacto con cierre de seguridad.

9) Peso máximo admitido 100Kg.

1) 4-wheel rollator with seat and basket that allows support for people with walking difficulties and thus be able to move as independently as possible.

2) Approved medical device based on the general requirements of safety, operation and provisions of Regulation (EU) 2017/745 of the European Parliament and the specific standard: UNE-EN ISO 11199-2:2021.

3) Hand brake with the possibility of locking rear wheels.

4) Locking big wheels with 3 positions 0-90-360°.

5) Height adjustable handle.

6) Removable bag with capacity for 25 litres.

7) Wide and ergonomic lateral support.

8) Compact folding with safety closure.

9) Maximum weight admitted 100 Kg.

Especificaciones Specifications

Colores Colours	211- Negro (Black) 739-Gris (Grey), 740-Granate (Maroon)
Material de la estructura Structure material	Acero, aluminio y plástico Steel, aluminium and plastic
Medidas de asiento [cm] Seat measures [cm]	34x35
Medidas exteriores - Ancho, Largo, Alto [cm] Outer dimensions - Width, Long, Height [cm]	52x63x91
Medidas exteriores plegado - Ancho, Largo, Alto [cm] Folded outer dimensions - Width, Long, Height [cm]	52x71x37
Tamaño de ruedas [cm] Wheels size [cm]	17
Medidas embalaje Ancho, Largo, Alto [cm] Packaging measures - Width, Long, Height [cm]	53x37x69
Peso [Kg] Weight [Kg]	7,5
Unidades embalaje Units per package	1
Peso embalaje bruto [Kg] Packaging gross weight [Kg]	9

