

Infrared thermometer

Model:YK-IRT2

Yonker is one of the biggest medical and healthcare product manufacturers in China.
 Range of product:Patient monitor, Pulse oximeter, Blood pressure monitor, Infrared Thermometer, Fetal Doppler,
 UV Radiation Treatment System, etc.The whole production line is medical CE certified and mostly FDA certified.

PRODUCT INTRODUCTION	FEATURES
Application It can measure human ear and forehead temperature, it suitable for medical unit and home use.	<ul style="list-style-type: none"> ● 1 second measure the temperature, easily and fast. ● Sensor measurement technology, high precision. ● Automatically power-off, if left idle for 60 seconds. ● One-key measurement, easy to use. ● Fever prompt, better to know your physical condition . ● Stores 34 sets recent measurement data, easy for your data contrast. ● Safety by infrared measuring, avoid the damage of the measuring by traditional mercury thermometer.
Contraindications None	
Measurement unit °C & °F can be switched by pressing a button	
Range 34.0°C~43.0°C	
Accuracy < 34.0°C and >43.0°C :±0.3°C, 34.0°C~ 43.0 °C :±0.2°C	
Measuring position Automatic forehead sensing	
User Group Age: Adult/Children -Fever patient -Other people who need to have their temperature taken -Doctor	
Normal operating conditions Ambient temperature: 15°C~40°C Relative humidity: 30%~85% Atmospheric pressure: 70kPa~106kPa	
Storage and transportation temperature Ambient temperature: -20°C~55°C Relative humidity: 15%~93% Atmospheric pressure: 50kPa~105kPa	
Type of battery DC 3V AAA alkaline battery *2	
Product dimension about 150mm×75mm×50mm (length × width × height)	
Product weight 120g (without batteries)	
Package including <ul style="list-style-type: none"> ● 1pcs host machine ● 1pcs manual for use ● 2pcs batteries 	



YK-IRT2

Infrared thermometer

MORE PRECISION MORE SECURITY

ENCYIRT2202307061ZZW

NOTE: The specifications are subject to changes without prior notice.

ACCURATE TEMPERATURE MEASUREMENT AND SAFETY CARE

YK-IRT2



01

Accuracy

Measuring sensor technology, high measurement accuracy

02

Intelligence

High temperature warning, sound and light warning

03

Intimate

Store data for easy analysis

04

Professional

Authoritative certification, quality assured



One second speed test, peace of mind is preferred

CLASSIC —Yonker Infrared thermometer YK-IRT2, accurate perception of each degree

Infrared temperature measurement, accurate capture

One second temperature measurement, accurate and fast. High-precision and high-sensitivity sensors with intelligent temperature algorithm can quickly capture the temperature, effectively reduce the measurement error, and operate more smoothly.

Big HD screen, quick reading

Large HD backlight screen, high brightness font, simple interface design, clear function display. It can also be read easily at night.

High temperature warning, sound and light warning

Intelligent warning, timely and intimate. Green, yellow and red backlight display indicates whether the body temperature is normal, and has a sound reminder, the measurement results are clear at a glance.

Data memory, easy comparison

One - key conversion mode, one - key check memory. Memory function, save the last 34 memory measurements (not memory Lo/Hi), real-time control of body temperature fluctuations, easy to analyze and compare.

Control the details, worry free

About 120g(net weight) body, light and portable, round appearance, comfortable to hold, authoritative certification, quality assurance.



EVERYWHERE, WE ALL DESIGN CAREFULLY.